

# RBcf: An VCF API for R.

Pierre Lindenbaum / @yokofakun/ Institut du Thorax . Nantes.

April 27, 2020

## 1 Abstract

RBcf uses the Htslib C API for parsing VCF and BCF files. This API was written by a regular user of the htsjdk library who doesn't like R.

A list of functions is available at: <https://github.com/lindenb/rbcf/blob/master/R/rbcf.R>

## 2 Examples

### 2.1 Show Htslib and Rbcf versions

Code:

```
# load the library
library(rbcf)
#print the version of the associated htslib
paste("HTSLIB:",htslib.version())
#print the version of rbcf
paste("RBCF:",rcbf.version())
```

Output:

```
[1] "HTSLIB: 1.10.2"
[1] "RBCF: 0.0-1"
```

### 2.2 Open and close a VCF file

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
```

```
fp <- bcf.open("./data/rotavirus_rf.01.vcf",FALSE)
# error (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# dispose the vcf reader
bcf.close(fp)
print("Done.")
```

Output:

```
[1] TRUE
[1] "Done."
```

## 2.3 Print the INFOs in the VCF header

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("./data/rotavirus_rf.01.vcf",FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# print INFO
bcf.infos(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
```

Output:

	ID	Number	Type
INDEL	INDEL	0	Flag
IDV	IDV	1	Integer
IMF	IMF	1	Float
DP	DP	1	Integer
VDB	VDB	1	Float
RPB	RPB	1	Float
MQB	MQB	1	Float
BQB	BQB	1	Float
MQSB	MQSB	1	Float
SGB	SGB	1	Float
MQOF	MQOF	1	Float
ICB	ICB	1	Float

HOB	HOB	1	Float	
AC	AC	A	Integer	
AN	AN	1	Integer	
DP4	DP4	4	Integer	
MQ	MQ	1	Integer	

  

INDEL		"Indicates that the v
IDV		"Maximum number of reads
IMF		"Maximum fraction of reads
DP		
VDB	"Variant Distance Bias for filtering splice-site artefacts in RNA-seq data	
RPB	"Mann-Whitney U test of Read Position Bias	
MQB	"Mann-Whitney U test of Mapping Quality Bias	
BQB	"Mann-Whitney U test of Base Quality Bias	
MQSB	"Mann-Whitney U test of Mapping Quality vs Strand Bias	
SGB	"Segreg	
MQOF	"Fraction of MQ0 reads	
ICB	"Inbreeding Coefficient Binomial test	
HOB	"Bias in the number of HOMs number	
AC	"Allele count in genotypes for each ALT allele, in the s	
AN	"Total number of alleles	
DP4	"Number of high-quality ref-forward , ref-reverse, alt-forward an	
MQ	"Aver	
[1]	TRUE	

## 2.4 Print the FORMATS in the VCF header

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("./data/rotavirus_rf.01.vcf",FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# print FORMAT
bcf.formats(fp)
# dispose the vcf reader
bcf.close(fp)
```

Output:



```

# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("./data/rotavirus_rf.01.vcf",FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# print the number of samples
paste("Number_of_samples:",bcf.nsamples(fp))
# get the name for the 1st sample
paste("First_sample:",bcf.sample.at(fp,1))
# get the 1-based index for the samples
bcf.sample2index(fp,c("S1","S2","S3","missing"))
# get all the samples
bcf.samples(fp)
# dispose the vcf reader
bcf.close(fp)

```

Output:

```

[1] "Number_of_samples:_5"
[1] "First_sample:_S1"
      S1      S2      S3 missing
      1      2      3      0
[1] "S1" "S2" "S3" "S4" "S5"
[1] TRUE

```

## 2.7 Print the Dictionary in the VCF header

Code:

```

# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("./data/rotavirus_rf.01.vcf",FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# print the dictionary
bcf.dictionary(fp)
# dispose the vcf reader
bcf.close(fp)

```

Output:

```

      chrom size
RF01  RF01 3302
RF02  RF02 2687
RF03  RF03 2592
RF04  RF04 2362
RF05  RF05 1579
RF06  RF06 1356
RF07  RF07 1074
RF08  RF08 1059
RF09  RF09 1062
RF10  RF10  751
RF11  RF11  666
[1] TRUE

```

## 2.8 Print the Indexed Chromosomes

Code:

```

# load rbcf
library(rbcf)
# Open the indexed VCF
fp <- bcf.open("./data/rotavirus_rf.02.vcf.gz")
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# get the indexed contigs
bcf.contigs(fp)
# dispose the vcf reader
bcf.close(fp)

```

Output:

```

[1] "RF01" "RF02" "RF03" "RF04" "RF05" "RF06" "RF07" "RF08" "RF09" "RF10"
[11] "RF11"
[1] TRUE

```

## 2.9 Scanning the variants

Code:

```

# load rbcf
library(rbcf)

```

```

# create a function counting variants in a VCF
count.variants<-function(filename) {
  # we don't need the index for this file
  fp <- bcf.open(filename,FALSE)
  # error on opening
  if(is.null(fp)) return(-1)
  # number of variants
  n<-0
  # loop while we can read a variant
  while(!is.null(vc<-bcf.next(fp))) {
    # increment the count
    n<-n+1
  }
  # dispose the vcf reader
  bcf.close(fp)
  # return the number of variant
  n
}

# filenames
vcfs<-c(
  "./data/gnomad.exomes.r2.0.1.sites.bcf",
  "./data/rotavirus_rf.01.vcf",
  "./data/rotavirus_rf.02.vcf.gz",
  "./data/rotavirus_rf.03.vcf.gz",
  "./data/rotavirus_rf.04.bcf"
)

# print the number of variants for each vcf
for(f in vcfs) {
  cat(paste(f," ",count.variants(f),"\n"))
}

```

Output:

```

./data/gnomad.exomes.r2.0.1.sites.bcf    50
./data/rotavirus_rf.01.vcf    45
./data/rotavirus_rf.02.vcf.gz    45
./data/rotavirus_rf.03.vcf.gz    45
./data/rotavirus_rf.04.bcf    45

```

## 2.10 Scanning the variants

Code:

```
# load rbcf
library(rbcf)

# create a function counting variants in a VCF
count.variants<-function(filename,predicate) {
  # we don't need the index for this file
  fp <- bcf.open(filename,FALSE)
  # error on opening
  if(is.null(fp)) return(-1)
  # number of variants
  n<-0
  # loop while we can read a variant
  while(!is.null(vc<-bcf.next(fp))) {
    # test the variant
    if(predicate(vc)) {
      # increment the count
      n<-n+1
    }
  }
  # dispose the vcf reader
  bcf.close(fp)
  # return the number of variant
  n
}

# A vcf
filename <- "../data/gnomad.exomes.r2.0.1.sites.bcf"
# filters
filters<-list(
  list("desc"="accept_all","predicate"=function(ctx) {TRUE} ),
  list("desc"="accept_none","predicate"=function(ctx) {FALSE} ),
  list("desc"="CHROM_is_1","predicate"=function(ctx) { variant.contig(ctx)%
  list("desc"="POS_is_even","predicate"=function(ctx) { (variant.pos(ctx)%
  list("desc"="PASS_filter","predicate"=function(ctx) {!variant.is.filtere
  list("desc"="count(FILTER)>1","predicate"=function(ctx) {length(variant.
  list("desc"="FILTER_contains_SEGDUP","predicate"=function(ctx) {variant.
  list("desc"="SNP","predicate"=function(ctx) {variant.is.snp(ctx)} ),
  list("desc"="POS!=END","predicate"=function(ctx) { variant.pos(ctx)!=var
```



```

list("desc"="not_diallelic","predicate"=function(ctx) {variant.nalleles(
list("desc"="REF_is_'A'", "predicate"=function(ctx) {variant.reference(ct
list("desc"="any_allele_is_'A'", "predicate"=function(ctx) {"A"  %in% var
list("desc"="any_ALT_allele_is_'A'", "predicate"=function(ctx) {"A"  %in%
list("desc"="No_QUAL", "predicate"=function(ctx) {!variant.has.qual(ctx)}
list("desc"="variant_has_ID", "predicate"=function(ctx) {variant.has.id(c
list("desc"="variant_ID_match_'rs1*'", "predicate"=function(ctx) {grepl(
list("desc"="variant_has_INFO/AF_NFE", "predicate"=function(ctx) {variant
list("desc"="variant_has_INFO/AF_NFE_>_1E-5", "predicate"=function(ctx) {
list("desc"="Missense_in_PLEKHN1_(VEP)", "predicate"=function(ctx) {
    # NO VEP annotation ?
    if(!variant.has.attribute(ctx,"CSQ")) return(FALSE);
    # get VEP annotation
    predictions <- variant.vep(ctx)
    # In SCN5A
    predictions <- predictions[which(predictions$SYMBOL=="PLEKHN1"),
    # Consequence must contain missense
    predictions <- predictions[grep("missense_variant",predictions$C
    nrow(predictions)>0
  })
)

# count the variant for each filter
for(flt in filters) {
  print(paste(basename(filename), "_filter:",flt[["desc"]], "_count:",count.
}

```

Output:

```

[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_accept_all_count:_50_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_accept_none_count:_0_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_CHROM_is_'1'_count:_50_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_POS_is_even_count:_24_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_PASS_filter_count:_48_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_count(FILTER)>1_count:_2_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_FILTER_contains_SEGDUP_count:_1_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_SNP_count:_47_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_POS!=END_count:_3_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_not_diallelic_count:_8_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_REF_is_'A'_count:_6_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_any_allele_is_'A'_count:_27_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_any_ALT_allele_is_'A'_count:_21_\n"

```

```
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_No_QUAL_count:_1_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_variant_has_ID_count:_34_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_variant_ID_match_'rs1*'_count:_2_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_variant_has_INFO/AF_NFE_count:_50_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_variant_has_INFO/AF_NFE>_1E-5_count:_1_\n"
[1] "gnomad.exomes.r2.0.1.sites.bcf_filter:_Missense_in_PLEKHN1_(VEP)_count:_1_\n"
```

## 2.11 Print a VEP table for a Variant

Code:

```
# load rbcf
library(rbcf)
# A vcf
filename <- "./data/gnomad.exomes.r2.0.1.sites.bcf"
# we don't need the index for this file
fp <- bcf.open(filename,FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# current variant
vc <- NULL
while(!is.null(vc<-bcf.next(fp))) {
  #find the first variant having an INFO/CSQ attribute
  if(variant.has.attribute(vc,"CSQ")) break;
}

if(!is.null(vc)) {
  # get the VEP table for the variant
  predictions<-variant.vep(vc)
}

# dispose the vcf reader
bcf.close(fp)
# show
predictions
```

Output:

```
[1] TRUE
      Allele      Consequence  IMPACT  SYMBOL  Gene
1      C downstream_gene_variant MODIFIER  KLHL17 ENSG00000187961
2      A downstream_gene_variant MODIFIER  KLHL17 ENSG00000187961
```

3	C	downstream_gene_variant	MODIFIER	C1orf170	ENSG00000187642
4	A	downstream_gene_variant	MODIFIER	C1orf170	ENSG00000187642
5	C	intron_variant	MODIFIER	PLEKHN1	ENSG00000187583
6	A	intron_variant	MODIFIER	PLEKHN1	ENSG00000187583
7	C	intron_variant	MODIFIER	PLEKHN1	ENSG00000187583
8	A	intron_variant	MODIFIER	PLEKHN1	ENSG00000187583
9	C	intron_variant	MODIFIER	PLEKHN1	ENSG00000187583
10	A	intron_variant	MODIFIER	PLEKHN1	ENSG00000187583
11	C	downstream_gene_variant	MODIFIER	C1orf170	ENSG00000187642
12	A	downstream_gene_variant	MODIFIER	C1orf170	ENSG00000187642
13	C	downstream_gene_variant	MODIFIER	C1orf170	ENSG00000187642
14	A	downstream_gene_variant	MODIFIER	C1orf170	ENSG00000187642
15	C	upstream_gene_variant	MODIFIER	PLEKHN1	ENSG00000187583
16	A	upstream_gene_variant	MODIFIER	PLEKHN1	ENSG00000187583
17	C	upstream_gene_variant	MODIFIER	PLEKHN1	ENSG00000187583
18	A	upstream_gene_variant	MODIFIER	PLEKHN1	ENSG00000187583

	Feature_type	Feature	BIOTYPE	EXON	INTRON
1	Transcript	ENST00000338591	protein_coding		
2	Transcript	ENST00000338591	protein_coding		
3	Transcript	ENST00000341290	protein_coding		
4	Transcript	ENST00000341290	protein_coding		
5	Transcript	ENST00000379407	protein_coding		2/14
6	Transcript	ENST00000379407	protein_coding		2/14
7	Transcript	ENST00000379409	protein_coding		2/14
8	Transcript	ENST00000379409	protein_coding		2/14
9	Transcript	ENST00000379410	protein_coding		2/15
10	Transcript	ENST00000379410	protein_coding		2/15
11	Transcript	ENST00000433179	protein_coding		
12	Transcript	ENST00000433179	protein_coding		
13	Transcript	ENST00000479361	retained_intron		
14	Transcript	ENST00000479361	retained_intron		
15	Transcript	ENST00000480267	retained_intron		
16	Transcript	ENST00000480267	retained_intron		
17	Transcript	ENST00000491024	protein_coding		
18	Transcript	ENST00000491024	protein_coding		

HGVSc HGVSp cDNA\_position CDS\_position

1  
2  
3  
4  
5 ENST00000379407.3:c.184-51G>C

6	ENST000000379407.3:c.184-51G>A								
7	ENST000000379409.2:c.184-51G>C								
8	ENST000000379409.2:c.184-51G>A								
9	ENST000000379410.3:c.184-51G>C								
10	ENST000000379410.3:c.184-51G>A								
11									
12									
13									
14									
15									
16									
17									
18									
	Protein_position	Amino_acids	Codons	Existing_variation	ALLELE_NUM	DISTANCE			
1				rs540662886	1	4511			
2				rs540662886	2	4511			
3				rs540662886	1	4978			
4				rs540662886	2	4978			
5				rs540662886	1				
6				rs540662886	2				
7				rs540662886	1				
8				rs540662886	2				
9				rs540662886	1				
10				rs540662886	2				
11				rs540662886	1	4973			
12				rs540662886	2	4973			
13				rs540662886	1	4979			
14				rs540662886	2	4979			
15				rs540662886	1	649			
16				rs540662886	2	649			
17				rs540662886	1	3286			
18				rs540662886	2	3286			
	STRAND	FLAGS	VARIANT_CLASS	MINIMISED	SYMBOL_SOURCE	HGNC_ID	CANONICAL		
1	1		SNV		HGNC	24023	YES		
2	1		SNV		HGNC	24023	YES		
3	-1		SNV		HGNC	28208			
4	-1		SNV		HGNC	28208			
5	1		SNV		HGNC	25284			
6	1		SNV		HGNC	25284			
7	1		SNV		HGNC	25284			
8	1		SNV		HGNC	25284			

9	1		SNV		HGNC	25284	YES	
10	1		SNV		HGNC	25284	YES	
11	-1		SNV		HGNC	28208	YES	
12	-1		SNV		HGNC	28208	YES	
13	-1		SNV		HGNC	28208		
14	-1		SNV		HGNC	28208		
15	1		SNV		HGNC	25284		
16	1		SNV		HGNC	25284		
17	1	cds_start_NF	SNV		HGNC	25284		
18	1	cds_start_NF	SNV		HGNC	25284		
	TSL	APPRIS	CCDS	ENSP	SWISSPROT	TREMBL	UNIPARC	
1			CCDS30550.1	ENSP00000343930	Q6TDP4	Q0VGE6&B3KXL7	UPI00001DFBF0	
2			CCDS30550.1	ENSP00000343930	Q6TDP4	Q0VGE6&B3KXL7	UPI00001DFBF0	
3				ENSP00000343864			UPI000022DAF4	
4				ENSP00000343864			UPI000022DAF4	
5			CCDS53256.1	ENSP00000368717	Q494U1	J3KSM5	UPI00005764FF	
6			CCDS53256.1	ENSP00000368717	Q494U1	J3KSM5	UPI00005764FF	
7				ENSP00000368719	Q494U1	J3KSM5	UPI0000D61E06	
8				ENSP00000368719	Q494U1	J3KSM5	UPI0000D61E06	
9			CCDS4.1	ENSP00000368720	Q494U1	J3KSM5	UPI00001416D8	
10			CCDS4.1	ENSP00000368720	Q494U1	J3KSM5	UPI00001416D8	
11				ENSP00000414022	Q5SV97		UPI0000418FB0	
12				ENSP00000414022	Q5SV97		UPI0000418FB0	
13								
14								
15								
16								
17				ENSP00000462558		J3KSM5	UPI0000268AE1F	
18				ENSP00000462558		J3KSM5	UPI0000268AE1F	
	GENE_PHENO	SIFT	PolyPhen	DOMAINS	HGVS_OFFSET	GMAF	AFR_MAF	AMR_MAF
1						C:0.0008	C:0	C:0
2						C:0.0008	C:0	C:0
3						C:0.0008	C:0	C:0
4						C:0.0008	C:0	C:0
5						C:0.0008	C:0	C:0
6						C:0.0008	C:0	C:0
7						C:0.0008	C:0	C:0
8						C:0.0008	C:0	C:0
9						C:0.0008	C:0	C:0
10						C:0.0008	C:0	C:0
11						C:0.0008	C:0	C:0

12						C:0.0008	C:0	C:0
13						C:0.0008	C:0	C:0
14						C:0.0008	C:0	C:0
15						C:0.0008	C:0	C:0
16						C:0.0008	C:0	C:0
17						C:0.0008	C:0	C:0
18						C:0.0008	C:0	C:0
	EAS_MAF	EUR_MAF	SAS_MAF	AA_MAF	EA_MAF	ExAC_MAF	ExAC_Adj_MAF	ExAC_AFR_MAF
1	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
2	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
3	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
4	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
5	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
6	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
7	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
8	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
9	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
10	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
11	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
12	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
13	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
14	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
15	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
16	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
17	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
18	C:0	C:0.004	C:0	C:0			C:0	C:2.146e-04
	ExAC_AMR_MAF	ExAC_EAS_MAF	ExAC_FIN_MAF	ExAC_NFE_MAF	ExAC_OTH_MAF			
1		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
2		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
3		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
4		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
5		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
6		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
7		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
8		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
9		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
10		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
11		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
12		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
13		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	
14		C:0	C:0.0002281	C:0.002986		C:0	C:1.606e-05	

15	C:0	C:0.0002281	C:0.002986	C:0	C:1.606e-05			
16	C:0	C:0.0002281	C:0.002986	C:0	C:1.606e-05			
17	C:0	C:0.0002281	C:0.002986	C:0	C:1.606e-05			
18	C:0	C:0.0002281	C:0.002986	C:0	C:1.606e-05			
	ExAC_SAS_MAF	CLIN_SIG	SOMATIC	PHENO	PUBMED	MOTIF_NAME	MOTIF_POS	HIGH_INF_POS
1	C:0							
2	C:0							
3	C:0							
4	C:0							
5	C:0							
6	C:0							
7	C:0							
8	C:0							
9	C:0							
10	C:0							
11	C:0							
12	C:0							
13	C:0							
14	C:0							
15	C:0							
16	C:0							
17	C:0							
18	C:0							
	MOTIF_SCORE_CHANGE	LoF	LoF_filter	LoF_flags	LoF_info"			
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								

## 2.12 Print a SNPEFF table for a Variant

Code:

```
# load rbcf
library(rbcf)
# A vcf
filename <- "./data/rotavirus_rf.ann.vcf.gz"
# we don't need the index for this file
fp <- bcf.open(filename,FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# current variant
vc <- NULL
while(!is.null(vc<-bcf.next(fp))) {
  #find the first variant having an INFO/ANN attribute
  if(variant.has.attribute(vc,"ANN")) break;
}
if(!is.null(vc)) {
  # get SNPEFF table
  predictions<-variant.snpeff(vc)
}
# dispose the vcf reader
bcf.close(fp)
# show
predictions
```

Output:

```
[1] TRUE
  Allele      Annotation Annotation_Impact  Gene_Name  Gene_ID
1      C missense_variant      MODERATE Gene_18_3284 Gene_18_3284
  Feature_Type Feature_ID Transcript_BioType Rank  HGVS.c      HGVS.p
1  transcript AAA47319.1    protein_coding  1/1 c.952A>C p.Lys318Gln
  cDNA.pos / cDNA.length CDS.pos / CDS.length AA.pos / AA.length Distance
1              952/3267          952/3267          318/1088
  ERRORS / WARNINGS / INFO'"
1
```



## 2.13 Query the indexed vcf using intervals

Code:

```
# load rbcf
library(rbcf)

# create a function counting variants in a VCF, in some intervals
count.variants<-function(filename,intervals) {
  # open the indexed VCF
  fp <- bcf.open(filename)
  # error on opening
  if(is.null(fp)) return(-1)
  # loop over the intervals
  for(interval in intervals) {
    # try query the interval
    if(bcf.query(fp,interval)) {
      # number of variants
      n<-0
      # loop while we can read a variant
      while(!is.null(vc<-bcf.next(fp))) {
        # increment the count
        n<-n+1
      }
      print(paste("Number_of_variants_in_",basename(filename),
        })
      # query failed
    else {
      print(paste("Cannot_query_",basename(filename),"/'",interval))
    }
  }
  # dispose the vcf reader
  bcf.close(fp)
}

some_intervals <-c("", "RF03", "RF03:2000-3000", "1:1-10000000", "chr1")
count.variants("./data/rotavirus_rf.02.vcf.gz",some_intervals)
count.variants("./data/1000G.ALL.2of4intersection.20100804.genotypes.bcf",some_i

# another way to query is set collect=TRUE to return a vector of variant
fp <- bcf.open("./data/rotavirus_rf.02.vcf.gz")
print(paste("Number_of_variants_using_collect:",length(bcf.query(fp,"RF03",colle
```

```
bcf.close(fp)
```

Output:

```
[1] "Cannot_query_rotavirus_rf.02.vcf.gz/''"
[1] "Number_of_variants_in_rotavirus_rf.02.vcf.gz/'RF03':8"
[1] "Number_of_variants_in_rotavirus_rf.02.vcf.gz/'RF03:2000-3000':4"
[1] "Cannot_query_rotavirus_rf.02.vcf.gz/'1:1-10000000'"
[1] "Cannot_query_rotavirus_rf.02.vcf.gz/'chr1'"
[1] TRUE
[1] "Cannot_query_1000G.ALL.2of4intersection.20100804.genotypes.bcf/''"
[1] "Cannot_query_1000G.ALL.2of4intersection.20100804.genotypes.bcf/'RF03'"
[1] "Cannot_query_1000G.ALL.2of4intersection.20100804.genotypes.bcf/'RF03:2000-3000'"
[1] "Number_of_variants_in_1000G.ALL.2of4intersection.20100804.genotypes.bcf/'1:1-10000000'"
[1] "Cannot_query_1000G.ALL.2of4intersection.20100804.genotypes.bcf/'chr1'"
[1] TRUE
[1] "Number_of_variants_using_collect:8"
[1] TRUE
```

## 2.14 Attribute in INFO

Code:

```
# load rbcf
library(rbcf)

# find given variant
find.variant<-function(fp,contig,pos) {
  if(!bcf.query(fp,paste(contig,":",pos,"-",pos,sep=""))) return(NULL)
  # loop while we can read a variant
  while(!is.null(vc<-bcf.next(fp))) {
    return(vc)
  }
  return(NULL)
}

filename<-"/data/gnomad.exomes.r2.0.1.sites.bcf"
# open the VCF with index
fp <- bcf.open(filename)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)

ctx <-find.variant(fp,"1",905608)
```

```

stopifnot(variant.has.attribute(ctx,"CSQ"))
print(paste("CSQ(no_split)",variant.string.attribute(ctx,"CSQ",split=FALSE)))
print(paste("CSQ(split)",variant.string.attribute(ctx,"CSQ")))
stopifnot(variant.has.attribute(ctx,"AN_POPMAX"))
print(paste("AN_POPMAX:",variant.int.attribute(ctx,"AN_POPMAX")))
stopifnot(variant.has.attribute(ctx,"AF_POPMAX"))
print(paste("AF_POPMAX:",variant.float.attribute(ctx,"AF_POPMAX")))
print(paste("flag:VQSR_NEGATIVE_TRAIN_SITE:",variant.flag.attribute(ctx,"VQSR_NE
# dispose the vcf reader
bcf.close(fp)

```

Output:

```

[1] "CSQ(no_split)TT|downstream_gene_variant|MODIFIER|KLHL17|ENSG00000187961|Tr
[1] "CSQ(split)TT|downstream_gene_variant|MODIFIER|KLHL17|ENSG00000187961|Tran
[2] "CSQ(split)TTA|downstream_gene_variant|MODIFIER|KLHL17|ENSG00000187961|Tran
[3] "CSQ(split)TT|downstream_gene_variant|MODIFIER|C1orf170|ENSG00000187642|Tr
[4] "CSQ(split)TTA|downstream_gene_variant|MODIFIER|C1orf170|ENSG00000187642|Tr
[5] "CSQ(split)TT|intron_variant|MODIFIER|PLEKHN1|ENSG00000187583|Transcript|E
[6] "CSQ(split)TTA|intron_variant|MODIFIER|PLEKHN1|ENSG00000187583|Transcript|E
[7] "CSQ(split)TT|intron_variant|MODIFIER|PLEKHN1|ENSG00000187583|Transcript|E
[8] "CSQ(split)TTA|intron_variant|MODIFIER|PLEKHN1|ENSG00000187583|Transcript|E
[9] "CSQ(split)TT|intron_variant|MODIFIER|PLEKHN1|ENSG00000187583|Transcript|E
[10] "CSQ(split)TTA|intron_variant|MODIFIER|PLEKHN1|ENSG00000187583|Transcript|E
[11] "CSQ(split)TT|downstream_gene_variant|MODIFIER|C1orf170|ENSG00000187642|Tr
[12] "CSQ(split)TTA|downstream_gene_variant|MODIFIER|C1orf170|ENSG00000187642|Tr
[13] "CSQ(split)TT|downstream_gene_variant|MODIFIER|C1orf170|ENSG00000187642|Tr
[14] "CSQ(split)TTA|downstream_gene_variant|MODIFIER|C1orf170|ENSG00000187642|Tr
[15] "CSQ(split)TT|upstream_gene_variant|MODIFIER|PLEKHN1|ENSG00000187583|Trans
[16] "CSQ(split)TTA|upstream_gene_variant|MODIFIER|PLEKHN1|ENSG00000187583|Trans
[17] "CSQ(split)TT|upstream_gene_variant|MODIFIER|PLEKHN1|ENSG00000187583|Trans
[18] "CSQ(split)TTA|upstream_gene_variant|MODIFIER|PLEKHN1|ENSG00000187583|Trans
[1] "AN_POPMAX: 106408" "AN_POPMAX: 106408"
[1] "AF_POPMAX: 1.87955993169453e-05" "AF_POPMAX: 9.39778965403093e-06"
[1] "flag:VQSR_NEGATIVE_TRAIN_SITE: FALSE"
[1] TRUE

```

## 2.15 Working with Genotypes

Code:

```
# load rbcf
```

```

library(rbcf)

# find given variant
find.variant<-function(fp,contig,pos) {
  if(!bcf.query(fp,paste(contig,":",pos,"-",pos,sep=""))) return(NULL)
  # loop while we can read a variant
  while(!is.null(vc<-bcf.next(fp))) {
    return(vc)
  }
  return(NULL)
}

filename<-"./data/1000G.ALL.2of4intersection.20100804.genotypes.bcf"
# open the VCF with index
fp <- bcf.open(filename)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)
# find a variant
ctx <-find.variant(fp,"1",10583)
print(paste("Number_of_genotypes_",variant.nsamples(ctx)))
# get 10-th genotype
gt<-variant.genotype(ctx,10)
print(paste("sample_",genotype.sample(gt)))
# get genotype by name
gt<-variant.genotype(ctx,"NA18997")
print(paste("sample_",genotype.sample(gt)))
print(paste("alleles_",genotype.alleles.idx0(gt)))
print(paste("genotype_ploidy_",genotype.ploidy(gt)))
print(paste("genotype_is_hom_ref_",genotype.homref(gt)))
print(paste("genotype_is_het_",genotype.het(gt)))
print(paste("genotype_is_het-non-ref_",genotype.hetnonref(gt)))
print(paste("genotype_is_phased_",genotype.phased(gt)))
print(paste("genotype_is_no_call_",genotype.nocall(gt)))
print(paste("genotype_FORMAT/OG_",genotype.string.attribute(gt,"OG")))
print(paste("genotype_FORMAT/GQ_",genotype.int.attribute(gt,"GQ")))# hum spec
print(paste("genotype_has_GQ_",genotype.has.gq(gt)))
print(paste("genotype_GQ_",genotype.gq(gt)))
print(paste("genotype_has_DP_",genotype.has.dp(gt)))
print(paste("genotype_DP_",genotype.int.attribute(gt,"DP")))
print(paste("genotype_DP_",genotype.dp(gt)))
print(paste("genotype_has_PL_",genotype.has.pl(gt)))
print(paste("genotype_PL_",genotype.pl(gt)))

```

```

print(paste("genotype_has_AD?", genotype.has.ad(gt)))
print(paste("genotype_AD", genotype.ad(gt)))

# dispose the vcf reader
bcf.close(fp)

```

**Output:**

```

[1] "Number_of_genotypes_629"
[1] "sample_HG00120"
[1] "sample_NA18997"
[1] "alleles_0" "alleles_1"
[1] "genotype_ploidy_2"
[1] "genotype_is_hom_ref_FALSE"
[1] "genotype_is_het_TRUE"
[1] "genotype_is_het-non-ref_FALSE"
[1] "genotype_is_phased_TRUE"
[1] "genotype_is_no_call_FALSE"
[1] "genotype_FORMAT/OG_1/1"
[1] "genotype_FORMAT/GQ_"
[1] "genotype_has_GQ_FALSE"
[1] "genotype_GQ_-1"
[1] "genotype_has_DP_TRUE"
[1] "genotype_DP_1"
[1] "genotype_DP_1"
[1] "genotype_has_PL_FALSE"
[1] "genotype_PL_"
[1] "genotype_has_AD_TRUE"
[1] "genotype_AD_4" "genotype_AD_1"
[1] TRUE

```

## 2.16 Writing variants to a new VCF/BCF file

**Code:**

```

# load rbcf
library(rbcf)
# vcf input filename
filenamein = "./data/rotavirus_rf.01.vcf"
# output vcf filename. "-" is standard output
filenameout = "-"

```

```

fp <- bcf.open(filenamein,FALSE)
# error on opening (exit 0 for tests)
if(is.null(fp)) quit(save="no",status=0,runLast=FALSE)

# create a new VCF writer using the header from 'fp'
out <- bcf.new.writer(fp,filenameout)
# error on opening (exit 0 for tests)
if(is.null(out)) quit(save="no",status=0,runLast=FALSE)

# loop while we can read a variant
while(!is.null(vc<-bcf.next(fp))) {
  # only write POS%10==0
  if(variant.pos(vc)%10==0) {
    # write variant
    bcf.write.variant(out,vc);
  }
}
# dispose the vcf reader
bcf.close(fp)
# dispose the vcf rwriter
bcf.close(out);

```

Output:

```

[1] TRUE
##fileformat=VCFv4.2
##FILTER=<ID=PASS,Description="All filters passed">
##samtoolsVersion=1.3.1+htslib-1.3.1
##samtoolsCommand=samtools mpileup -Ou -f rotavirus_rf.fa S1.bam S2.bam S3.bam
##reference=file://rotavirus_rf.fa
##contig=<ID=RF01,length=3302>
##contig=<ID=RF02,length=2687>
##contig=<ID=RF03,length=2592>
##contig=<ID=RF04,length=2362>
##contig=<ID=RF05,length=1579>
##contig=<ID=RF06,length=1356>
##contig=<ID=RF07,length=1074>
##contig=<ID=RF08,length=1059>
##contig=<ID=RF09,length=1062>
##contig=<ID=RF10,length=751>
##contig=<ID=RF11,length=666>
##ALT=<ID=*,Description="Represents allele(s) other than observed.">

```

```

##INFO=<ID=INDEL,Number=0,Type=Flag,Description="Indicates that the variant is
##INFO=<ID=IDV,Number=1,Type=Integer,Description="Maximum number of reads supp
##INFO=<ID=IMF,Number=1,Type=Float,Description="Maximum fraction of reads supp
##INFO=<ID=DP,Number=1,Type=Integer,Description="Raw read depth">
##INFO=<ID=VDB,Number=1,Type=Float,Description="Variant Distance Bias for filt
##INFO=<ID=RPB,Number=1,Type=Float,Description="Mann-Whitney U test of Read Pos
##INFO=<ID=MQB,Number=1,Type=Float,Description="Mann-Whitney U test of Mapping
##INFO=<ID=BQB,Number=1,Type=Float,Description="Mann-Whitney U test of Base Qu
##INFO=<ID=MQSB,Number=1,Type=Float,Description="Mann-Whitney U test of Mapping
##INFO=<ID=SGB,Number=1,Type=Float,Description="Segregation based metric.">
##INFO=<ID=MQOF,Number=1,Type=Float,Description="Fraction of MQ0 reads (smaller
##FORMAT=<ID=PL,Number=G,Type=Integer,Description="List of Phred-scaled genotyp
##FORMAT=<ID=GT,Number=1,Type=String,Description="Genotype">
##INFO=<ID=ICB,Number=1,Type=Float,Description="Inbreeding Coefficient Binomial
##INFO=<ID=HOB,Number=1,Type=Float,Description="Bias in the number of HOMs num
##INFO=<ID=AC,Number=A,Type=Integer,Description="Allele count in genotypes for
##INFO=<ID=AN,Number=1,Type=Integer,Description="Total number of alleles in ca
##INFO=<ID=DP4,Number=4,Type=Integer,Description="Number of high-quality ref-f
##INFO=<ID=MQ,Number=1,Type=Integer,Description="Average mapping quality">
##bcftools_callVersion=1.3-10-g820e1d6+htslib-1.2.1-267-g87141ea
##bcftools_callCommand=call -um -Oz -o rotavirus_rf.vcf.gz -
##bcftools_viewVersion=1.10-6-g2782d9f+htslib-1.2.1-1336-g7c16b56-dirty
##bcftools_viewCommand=view /home/lindenb/src/jvarkit/src/test/resources/rotavir
#CHROM POS ID REF ALT QUAL FILTER INFO FORMAT S1
RF01 970 . A C 48.6696 . DP=36;VDB=0.693968;SGB=1
RF03 2150 . T A 6.90687 . DP=37;VDB=0.557348;SGB=-
RF04 1900 . A C 36.8224 . DP=39;VDB=0.706942;SGB=7
RF04 1920 . A T 42.014 . DP=39;VDB=0.966939;SGB=0

```