

# RBcf: An VCF API for R.

Pierre Lindenbaum / yokofakun/ Institut du Thorax . Nantes.

April 23, 2020

## 1 Abstract

RBcf uses the Htslib C API for parsing VCF and BCF files. This API was written by a regular user of the htsjdk library who doesn't like R.

## 2 Examples

### 2.1 Open and close a VCF file

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
# dispose the vcf reader
bcf.close(fp)
```

Output:

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
>
```

## 2.2 Print the INFOs in the VCF header

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
info <- bcf.infos(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
info
```

Output:

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
> info <- bcf.infos(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> info
```

	ID	Number	Type
INDEL	INDEL	0	Flag
IDV	IDV	1	Integer
IMF	IMF	1	Float
DP	DP	1	Integer
VDB	VDB	1	Float
RPB	RPB	1	Float
MQB	MQB	1	Float
BQB	BQB	1	Float
MQSB	MQSB	1	Float
SGB	SGB	1	Float
MQOF	MQOF	1	Float
ICB	ICB	1	Float
HOB	HOB	1	Float
AC	AC	A	Integer
AN	AN	1	Integer

DP4	DP4	4 Integer
MQ	MQ	1 Integer

  

INDEL	"Indicates that the v
IDV	"Maximum number of reads
IMF	"Maximum fraction of reads
DP	
VDB	"Variant Distance Bias for filtering splice-site artefacts in RNA-seq data
RPB	"Mann-Whitney U test of Read Position Bias
MQB	"Mann-Whitney U test of Mapping Quality Bias
BQB	"Mann-Whitney U test of Base Quality Bias
MQSB	"Mann-Whitney U test of Mapping Quality vs Strand Bias
SGB	"Segreg
MQOF	"Fraction of MQ0 reads
ICB	"Inbreeding Coefficient Binomial test
HOB	"Bias in the number of HOMs number
AC	"Allele count in genotypes for each ALT allele, in the s
AN	"Total number of alleles
DP4	"Number of high-quality ref-forward , ref-reverse, alt-forward an
MQ	"Aver
>	

## 2.3 Print the FORMATS in the VCF header

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
fmts <- bcf.formats(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
fmts
```

Output:

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
```

```

> fmts <- bcf.formats(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> fmts
      ID Number      Type      Description
PL PL          G Integer "List of Phred-scaled genotype likelihoods"
GT GT          1 String   "Genotype"
>

```

## 2.4 Print the FILTERs in the VCF header

Code:

```

# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/gnomad.exomes.r2.0.1.sites.vcf",FALSE)
flt <- bcf.filters(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
flt

```

Output:

```

> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/gnomad.exomes.r2.0.1.sites.vcf",FALSE)
> flt <- bcf.filters(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> flt

```

	ID
PASS	PASS
AC0	AC0
InbreedingCoeff	InbreedingCoeff
LCR	LCR

```

RF
SEGDUPE
RF
SEGDUPE

PASS
AC0 "Allele Count is zero (i.e. no high-confidence genotype (GQ >= 2
InbreedingCoeff
LCR
RF "Failed random forests filter
SEGDUPE
>

```

## 2.5 Print the Samples in the VCF header

Code:

```

# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
# print the number of samples
cat(paste("Num. Samples=",bcf.nsamples(fp),".\n"))
# get the name for the 1st sample
cat(paste("First sample is ",bcf.sample1(fp,1),".\n"))

# get the samples
samples <- bcf.samples(fp)
# dispose the vcf reader
bcf.close(fp)
# print the list
samples

```

Output:

```

> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
> # print the number of samples
> cat(paste("Num. Samples=",bcf.nsamples(fp),".\n"))
Num. Samples= 5 .
> # get the name for the 1st sample
> cat(paste("First sample is ",bcf.sample1(fp,1),".\n"))

```

```

First sample is S1 .
>
> # get the samples
> samples <- bcf.samples(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the list
> samples
[1] "S1" "S2" "S3" "S4" "S5"
>

```

## 2.6 Print the Dictionary in the VCF header

Code:

```

# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
dict <- bcf.dictionary(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
dict

```

Output:

```

> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
> dict <- bcf.dictionary(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> dict
      chrom size
RF01  RF01 3302
RF02  RF02 2687
RF03  RF03 2592

```

```

RF04 RF04 2362
RF05 RF05 1579
RF06 RF06 1356
RF07 RF07 1074
RF08 RF08 1059
RF09 RF09 1062
RF10 RF10 751
RF11 RF11 666
>

```

## 2.7 Print the Indexed Chromosomes

Code:

```

# load rbcf
library(rbcf)
# Open the indexed VCF
fp <- bcf.open("../tests/data/rotavirus_rf.02.vcf.gz")
# get the indexed contigs
contigs <- bcf.contigs(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
contigs

```

Output:

```

> # load rbcf
> library(rbcf)
> # Open the indexed VCF
> fp <- bcf.open("../tests/data/rotavirus_rf.02.vcf.gz")
> # get the indexed contigs
> contigs <- bcf.contigs(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> contigs
[1] "RF01" "RF02" "RF03" "RF04" "RF05" "RF06" "RF07" "RF08" "RF09" "RF10"
[11] "RF11"
>

```