RBcf: An VCF API for R.

Pierre Lindenbaum / @yokofakun/ Institut du Thorax . Nantes.

April 24, 2020

1 Abstract

RBcf uses the Htslib C API for parsing VCF and BCF files. This API was written by a regular user of the htsjdk library who doesn't like R.

A list of functions is available at: https://github.com/lindenb/rbcf/blob/master/R/rbcf.R

2 Examples

2.1 Htslib and Rbcf versions

Code:

```
# load the library
library(rbcf)
#print the version of the associated htslib
htslib.version()
#print the version of rbcf
rcbf.version()
```

```
> # load the library
> library(rbcf)
> #print the version of the associated htslib
> htslib.version()
[1] "1.10.2-dirty"
> #print the version of rbcf
> rcbf.version()
[1] "0.0-1"
>
```

2.2 Open and close a VCF file

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
# dispose the vcf reader
bcf.close(fp)</pre>
```

Output:

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
```

2.3 Print the INFOs in the VCF header

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
info <- bcf.infos(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
info</pre>
```

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
> info <- bcf.infos(fp)
> # dispose the vcf reader
> bcf.close(fp)
```

```
[1] TRUE
> # print the table
> info
                               ID Number
                                                                            Type
INDEL INDEL
                                                          0
                                                                            Flag
IDV
                           IDV
                                                          1 Integer
IMF
                           IMF
                                                          1
                                                                        Float
DP
                               DP
                                                          1 Integer
VDB
                           VDB
                                                          1
                                                                        Float
RPB
                           RPB
                                                          1
                                                                        Float
MQB
                           MQB
                                                                       Float
                                                          1
BQB
                           BQB
                                                          1
                                                                      Float
MQSB
                        MQSB
                                                          1
                                                                      Float
SGB
                           SGB
                                                          1
                                                                      Float
MQOF
                        MQOF
                                                          1
                                                                       Float
ICB
                           ICB
                                                          1
                                                                       Float
HOB
                           HOB
                                                                        Float
                                                          1
AC
                              AC
                                                          A Integer
AN
                              AN
                                                          1 Integer
DP4
                           DP4
                                                          4 Integer
MQ
                              MQ
                                                          1 Integer
INDEL
                                                                                                                                                                                                           "Indicates_that_the_v
IDV
                                                                                                                                                                                             "Maximum_number_of_reads_
IMF
                                                                                                                                                                                      "Maximum<sub>□</sub>fraction<sub>□</sub>of<sub>□</sub>reads<sub>□</sub>
DP
VDB
                     "Variant_Distance_Bias_for_filtering_splice-site_artefacts_in_RNA-seq_data
RPB
                                                                                                                                   "Mann-Whitney Uutest of Read Position Bias
MQB
                                                                                                                            "Mann-Whitney Uutest of Mapping Quality Bias
BQB
                                                                                                                                      "Mann-Whitney Uutest of Base Quality Bias
MQSB
                                                                                         \verb|"Mann-Whitney|| U_{\sqcup} test_{\sqcup} of_{\sqcup} \verb|Mapping|| Quality_{\sqcup} vs_{\sqcup} Strand_{\sqcup} Bias
SGB
                                                                                                                                                                                                                                                           "Segreg
MQOF
                                                                                                                                                                                                    "Fraction of MQ0 reads to make the reads of the large term of the
ICB
                                                                                                                                                    "Inbreeding Coefficient Binomial test
HOB
                                                                                                                                                           "Bias_in_the_number_of_HOMs_number_
AC
                                                                                   "Allele\_count\_in\_genotypes\_for\_each\_ALT\_allele,\_in\_the\_s
                                                                                                                                                                                             "Total_{\sqcup}number_{\sqcup}of_{\sqcup}alleles_{\sqcup}
AN
DP4
                                                    \verb|"Number_lof_lhigh-quality_ref-forward_l,_lref-reverse,_lalt-forward_lan| \\
MQ
                                                                                                                                                                                                                                                                 "Aver
```

2.4 Print the FORMATs in the VCF header

Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
fmts <- bcf.formats(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
fmts</pre>
```

Output:

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)</pre>
> fmts <- bcf.formats(fp)</pre>
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> fmts
   ID Number
                Type
                                                        Description
PL PL
           G Integer "List⊔of⊔Phred-scaled⊔genotype⊔likelihoods"
           1 String
GT GT
                                                         "Genotype"
```

2.5 Print the FILTERs in the VCF header

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/gnomad.exomes.r2.0.1.sites.vcf",FALSE)
flt <- bcf.filters(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
flt</pre>
```

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/gnomad.exomes.r2.0.1.sites.vcf",FALSE)</pre>
> flt <- bcf.filters(fp)</pre>
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> flt
                                ID
PASS
                              PASS
AC0
                               ACO
InbreedingCoeff InbreedingCoeff
LCR
                               LCR
RF
                                RF
SEGDUP
                           SEGDUP
PASS
ACO
                 "Allele_Count_is_zero_(i.e._no_high-confidence_genotype_(GQ_>=_2
InbreedingCoeff
LCR
                                                          "Failed_{\sqcup}random_{\sqcup}forests_{\sqcup}filte
RF
SEGDUP
```

2.6 Print the Samples in the VCF header

The samples are defined in the '#CHROM' line of the VCF Code:

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
# print the number of samples
bcf.nsamples(fp)
# get the name for the 1st sample
bcf.sample.at(fp,1)
# get the 1-based index for the samples
bcf.sample2index(fp,c("S1","S2","S3","missing"))</pre>
```

```
# get all the samples
bcf.samples(fp)
# dispose the vcf reader
bcf.close(fp)
```

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)</pre>
> # print the number of samples
> bcf.nsamples(fp)
[1] 5
> # get the name for the 1st sample
> bcf.sample.at(fp,1)
[1] "S1"
> # get the 1-based index for the samples
> bcf.sample2index(fp,c("S1","S2","S3","missing"))
    S1
            S2 S3 missing
     1
             2
                     3
                              0
> # get all the samples
> bcf.samples(fp)
[1] "S1" "S2" "S3" "S4" "S5"
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
```

2.7 Print the Dictionary in the VCF header

```
# load rbcf
library(rbcf)
# we don't need the index for this file
fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)
dict <- bcf.dictionary(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
dict</pre>
```

```
> # load rbcf
> library(rbcf)
> # we don't need the index for this file
> fp <- bcf.open("../tests/data/rotavirus_rf.01.vcf",FALSE)</pre>
> dict <- bcf.dictionary(fp)</pre>
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> dict
     chrom size
RF01 RF01 3302
RF02 RF02 2687
RF03 RF03 2592
RF04 RF04 2362
RF05 RF05 1579
RF06 RF06 1356
RF07 RF07 1074
RF08 RF08 1059
RF09 RF09 1062
RF10 RF10 751
RF11 RF11 666
```

2.8 Print the Indexed Chromosomes

Code:

```
# load rbcf
library(rbcf)
# Open the indexed VCF
fp <- bcf.open("../tests/data/rotavirus_rf.02.vcf.gz")
# get the indexed contigs
contigs <- bcf.contigs(fp)
# dispose the vcf reader
bcf.close(fp)
# print the table
contigs</pre>
```

```
> # load rbcf
> library(rbcf)
> # Open the indexed VCF
> fp <- bcf.open("../tests/data/rotavirus_rf.02.vcf.gz")
> # get the indexed contigs
> contigs <- bcf.contigs(fp)
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # print the table
> contigs
[1] "RF01" "RF02" "RF03" "RF04" "RF05" "RF06" "RF07" "RF08" "RF09" "RF10"
[11] "RF11"
>
```

2.9 Scanning the variants

```
# load rbcf
library(rbcf)
# create a function counting variants in a VCF
count.variants<-function(filename) {</pre>
        # we don't need the index for this file
        fp <- bcf.open(filename,FALSE)</pre>
        # number of variants
        n<-0
        # loop while we can read a variant
        while(!is.null(vc<-bcf.next(fp))) {</pre>
                # increment the count
                n<-n+1
        # dispose the vcf reader
        bcf.close(fp)
        # return the number of variant
        n
# filenames
vcfs<-c(
```

```
"../tests/data/gnomad.exomes.r2.0.1.sites.vcf",
    "../tests/data/rotavirus_rf.01.vcf",
    "../tests/data/rotavirus_rf.02.vcf.gz",
    "../tests/data/rotavirus_rf.03.vcf.gz",
    "../tests/data/rotavirus_rf.04.bcf"
    )
# print the number of variants for each vcf
for(f in vcfs) {
    cat(paste(f,"\",count.variants(f),"\n"))
    }
```

```
> # load rbcf
> library(rbcf)
> # create a function counting variants in a VCF
> count.variants<-function(filename) {</pre>
        # we don't need the index for this file
        fp <- bcf.open(filename,FALSE)</pre>
        # number of variants
        n < -0
        # loop while we can read a variant
        while(!is.null(vc<-bcf.next(fp))) {</pre>
                # increment the count
                n<-n+1
        }
        # dispose the vcf reader
        bcf.close(fp)
        # return the number of variant
        n
+ }
> # filenames
> vcfs<-c(</pre>
        "../tests/data/gnomad.exomes.r2.0.1.sites.vcf",
        "../tests/data/rotavirus_rf.01.vcf",
        "../tests/data/rotavirus_rf.02.vcf.gz",
        "../tests/data/rotavirus_rf.03.vcf.gz",
        "../tests/data/rotavirus_rf.04.bcf"
> # print the number of variants for each vcf
```

```
> for(f in vcfs) {
+     cat(paste(f,"u",count.variants(f),"\n"))
+     }
../tests/data/gnomad.exomes.r2.0.1.sites.vcf 50
../tests/data/rotavirus_rf.01.vcf 45
../tests/data/rotavirus_rf.02.vcf.gz 45
../tests/data/rotavirus_rf.03.vcf.gz 45
../tests/data/rotavirus_rf.04.bcf 45
>
```

2.10 Scanning the variants

```
# load rbcf
library(rbcf)
# create a function counting variants in a VCF
count.variants<-function(filename,predicate) {</pre>
        # we don't need the index for this file
        fp <- bcf.open(filename,FALSE)</pre>
        # number of variants
        n<-0
        # loop while we can read a variant
        while(!is.null(vc<-bcf.next(fp))) {</pre>
                 # test the variant
                if(predicate(vc)) {
                         # increment the count
                         n < -n + 1
                         }
        }
        # dispose the vcf reader
        bcf.close(fp)
        # return the number of variant
# A vcf
filename <- "../tests/data/gnomad.exomes.r2.0.1.sites.vcf"
# filters
filters<-list(
```

```
list("desc"="accept_all", "predicate"=function(ctx) {TRUE} ),
        list("desc"="accept_none", "predicate"=function(ctx) {FALSE} ),
        list("desc"="CHROM_is_'1', "predicate"=function(ctx) { variant.contig(ct
        list("desc"="POS<sub>□</sub>is<sub>□</sub>even", "predicate"=function(ctx) { (variant.pos(ctx)%
        list("desc"="PASS_filter", "predicate"=function(ctx) {!variant.is.filtere
        list("desc"="count(FILTER)>1","predicate"=function(ctx) {length(variant.
        list("desc"="FILTER_contains_SEGDUP", "predicate"=function(ctx) {variant.
        list("desc"="SNP", "predicate"=function(ctx) {variant.is.snp(ctx)} ),
        list("desc"="POS!=END", "predicate"=function(ctx) { variant.pos(ctx)!=var
        list("desc"="not_diallelic", "predicate"=function(ctx) {variant.nalleles(
        list("desc"="REF<sub>□</sub>is<sub>□</sub>'A'", "predicate"=function(ctx) {variant.reference(ct
        list("desc"="anyualleleuisu'A'", "predicate"=function(ctx) {"A" %in% var
        list("desc"="any LALT allele is 'A'", "predicate"=function(ctx) {"A" %in%
        list("desc"="No⊔QUAL", "predicate"=function(ctx) {!variant.has.qual(ctx)}
        list("desc"="variant_has_ID", "predicate"=function(ctx) {variant.has.id(c
        list("desc"="variant_ID_match_'rs1*'_", "predicate"=function(ctx) {grepl(
        list("desc"="variant_has_INFO/AF_NFE", "predicate"=function(ctx) {variant
        list("desc"="variant_has_INFO/AF_NFE_>_1E-5", "predicate"=function(ctx) {
        list("desc"="Missense_{\sqcup}in_{\sqcup}PLEKHN1_{\sqcup}(VEP)", "predicate"=function(ctx) \ \{ list("desc"="Missense_{\sqcup}in_{\sqcup}PLEKHN1_{\sqcup}(VEP)", "predicate"=function(ctx) \ \} 
                 # NO VEP annotation ?
                 if(!variant.has.attribute(ctx, "CSQ")) return(FALSE);
                  # get VEP annotation
                 predictions <- variant.vep(ctx)</pre>
                 # In SCN5A
                 predictions <- predictions[which(predictions$SYMBOL=="PLEKHN1"),</pre>
                 # Consequence must contain missense
                 predictions <- predictions[grep("missense_variant",predictions$C</pre>
                 nrow(predictions)>0
                 })
         )
# count the variant for each filter
for(flt in filters) {
         cat(paste(flt[["desc"]],"\",count.variants(filename,flt[["predicate"]]),
Output:
> # load rbcf
> library(rbcf)
```

> # create a function counting variants in a VCF

```
count.variants<-function(filename, predicate) {</pre>
        # we don't need the index for this file
        fp <- bcf.open(filename,FALSE)</pre>
        # number of variants
        n<-0
        # loop while we can read a variant
        while(!is.null(vc<-bcf.next(fp))) {</pre>
                 # test the variant
                 if(predicate(vc)) {
                         # increment the count
                         n<-n+1
                         }
        }
        # dispose the vcf reader
        bcf.close(fp)
        # return the number of variant
+ }
> # A vcf
> filename <- "../tests/data/gnomad.exomes.r2.0.1.sites.vcf"</pre>
> # filters
> filters<-list(
        list("desc"="accept_all", "predicate"=function(ctx) {TRUE} ),
        list("desc"="accept_none", "predicate"=function(ctx) {FALSE} ),
        list("desc"="CHROM_is_''1', "predicate"=function(ctx) { variant.contig(ct
        list("desc"="POS<sub>□</sub>is<sub>□</sub>even", "predicate"=function(ctx) { (variant.pos(ctx)%
        list("desc"="PASS_filter", "predicate"=function(ctx) {!variant.is.filtere
        list("desc"="count(FILTER)>1", "predicate"=function(ctx) {length(variant.
        list("desc"="FILTER_{\sqcup}contains_{\sqcup}SEGDUP", "predicate"=function(ctx) \ \{variant.\}
        list("desc"="SNP","predicate"=function(ctx) {variant.is.snp(ctx)} ),
        list("desc"="POS!=END", "predicate"=function(ctx) { variant.pos(ctx)!=var
        list("desc"="not_diallelic", "predicate"=function(ctx) {variant.nalleles(
        list("desc"="REF<sub>□</sub>is<sub>□</sub>'A'", "predicate"=function(ctx) {variant.reference(ct
        list("desc"="anyualleleuisu'A'", "predicate"=function(ctx) {"A" %in% var
        list("desc"="any ALT allele is, 'A'", "predicate"=function(ctx) {"A" %in%
        list("desc"="No_QUAL", "predicate"=function(ctx) {!variant.has.qual(ctx)}
        list("desc"="variant_has_ID", "predicate"=function(ctx) {variant.has.id(c
        list("desc"="variant_|ID_match_'rs1*'_", "predicate"=function(ctx) {grepl(
        list("desc"="variant_has_INFO/AF_NFE", "predicate"=function(ctx) {variant
        list("desc"="variant_has_INFO/AF_NFE_>_1E-5", "predicate"=function(ctx) {
```

```
list("desc"="Missense_in_PLEKHN1_(VEP)", "predicate"=function(ctx) {
                # NO VEP annotation ?
                if(!variant.has.attribute(ctx, "CSQ")) return(FALSE);
                # get VEP annotation
                predictions <- variant.vep(ctx)</pre>
                # In SCN5A
                predictions <- predictions[which(predictions$SYMBOL=="PLEKHN1"),</pre>
                # Consequence must contain missense
                predictions <- predictions[grep("missense_variant",predictions$C</pre>
                nrow(predictions)>0
                })
        )
> # count the variant for each filter
> for(flt in filters) {
        cat(paste(flt[["desc"]],"_\",count.variants(filename,flt[["predicate"]]),
accept all
             50
accept none
CHROM is '1'
              50
POS is even
              24
PASS filter
count(FILTER)>1
FILTER contains SEGDUP
SNP 47
POS!=END
not diallelic
REF is 'A' 6
any allele is 'A'
                    27
any ALT allele is 'A'
                        21
No QUAL
         1
variant has ID
                 34
variant ID match 'rs1*'
variant has INFO/AF_NFE
                           50
variant has INFO/AF_NFE > 1E-5
                                  14
Missense in PLEKHN1 (VEP)
```

2.11 Print a VEP table for a Variant

Code:

```
# load rbcf
library(rbcf)
# A vcf
filename <- "../tests/data/gnomad.exomes.r2.0.1.sites.vcf"
# we don't need the index for this file
fp <- bcf.open(filename,FALSE)</pre>
# current variant
vc <- NULL
while(!is.null(vc<-bcf.next(fp))) {</pre>
        #find the first variant having an INFO/CSQ attribute
        if(variant.has.attribute(vc, "CSQ")) break;
        }
if(!is.null(vc)) {
        # get the VEP table for the variant
        predictions<-variant.vep(vc)</pre>
# dispose the vcf reader
bcf.close(fp)
# show
predictions
```

```
> # load rbcf
> library(rbcf)
> # A vcf
> filename <- "../tests/data/gnomad.exomes.r2.0.1.sites.vcf"
> # we don't need the index for this file
> fp <- bcf.open(filename, FALSE)
> # current variant
> vc <- NULL
> while(!is.null(vc<-bcf.next(fp))) {
+     #find the first variant having an INFO/CSQ attribute
+     if(variant.has.attribute(vc, "CSQ")) break;
+     }
> if(!is.null(vc)) {
```

```
# get the VEP table for the variant
        predictions<-variant.vep(vc)</pre>
        }
 # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # show
> predictions
   Allele
                      Consequence
                                     IMPACT
                                              SYMBOL
                                                                 Gene
        C downstream_gene_variant MODIFIER
                                              KLHL17 ENSG00000187961
2
        A downstream_gene_variant MODIFIER
                                              KLHL17 ENSG00000187961
3
        C downstream_gene_variant MODIFIER C1orf170 ENSG00000187642
4
        A downstream_gene_variant MODIFIER Clorf170 ENSG00000187642
5
                   intron variant MODIFIER PLEKHN1 ENSG00000187583
        C
6
        Α
                   intron_variant MODIFIER PLEKHN1 ENSG00000187583
7
                   intron_variant MODIFIER PLEKHN1 ENSG00000187583
        C
8
                   intron_variant MODIFIER PLEKHN1 ENSG00000187583
        Α
9
                   intron_variant MODIFIER PLEKHN1 ENSG00000187583
        C
                   intron_variant MODIFIER PLEKHN1 ENSG00000187583
10
        Α
11
        C downstream_gene_variant MODIFIER C1orf170 ENSG00000187642
        A downstream_gene_variant MODIFIER Clorf170 ENSG00000187642
12
13
        C downstream_gene_variant MODIFIER Clorf170 ENSG00000187642
        A downstream_gene_variant MODIFIER Clorf170 ENSG00000187642
14
15
            upstream_gene_variant MODIFIER PLEKHN1 ENSG00000187583
16
        Α
            upstream_gene_variant MODIFIER PLEKHN1 ENSG00000187583
17
            upstream_gene_variant MODIFIER
                                           PLEKHN1 ENSG00000187583
        C
18
            upstream_gene_variant MODIFIER PLEKHN1 ENSG00000187583
                        Feature
                                         BIOTYPE EXON INTRON
   Feature_type
1
     Transcript ENST00000338591
                                 protein_coding
2
     Transcript ENST00000338591
                                 protein_coding
3
                                 protein_coding
     Transcript ENST00000341290
4
     Transcript ENST00000341290
                                 protein_coding
5
     Transcript ENST00000379407
                                  protein_coding
                                                        2/14
6
     Transcript ENST00000379407
                                  protein_coding
                                                        2/14
     Transcript ENST00000379409
                                 protein_coding
                                                        2/14
8
     Transcript ENST00000379409
                                  protein_coding
                                                        2/14
9
     Transcript ENST00000379410
                                 protein_coding
                                                        2/15
10
     Transcript ENST00000379410
                                  protein_coding
                                                        2/15
11
     Transcript ENST00000433179
                                  protein_coding
12
     Transcript ENST00000433179
                                 protein_coding
```

```
13
     Transcript ENST00000479361 retained_intron
14
     Transcript ENST00000479361 retained_intron
15
     Transcript ENST00000480267 retained_intron
16
     Transcript ENST00000480267 retained_intron
17
     Transcript ENST00000491024
                                 protein_coding
18
     Transcript ENST00000491024 protein_coding
                            HGVSc HGVSp cDNA_position CDS_position
2
3
4
5
  ENST00000379407.3:c.184-51G>C
  ENST00000379407.3:c.184-51G>A
  ENST00000379409.2:c.184-51G>C
 ENST00000379409.2:c.184-51G>A
9 ENST00000379410.3:c.184-51G>C
10 ENST00000379410.3:c.184-51G>A
11
12
13
14
15
16
17
18
   Protein_position Amino_acids Codons Existing_variation ALLELE_NUM DISTANCE
                                                rs540662886
1
                                                                      1
                                                                            4511
2
                                                                      2
                                                rs540662886
                                                                            4511
3
                                                rs540662886
                                                                      1
                                                                            4978
4
                                                rs540662886
                                                                      2
                                                                            4978
5
                                                rs540662886
                                                                      1
6
                                                rs540662886
                                                                      2
7
                                                rs540662886
                                                                      1
8
                                                                      2
                                                rs540662886
9
                                                rs540662886
                                                                      1
10
                                                                      2
                                                rs540662886
11
                                                rs540662886
                                                                            4973
                                                                      1
12
                                                rs540662886
                                                                      2
                                                                            4973
13
                                                rs540662886
                                                                      1
                                                                            4979
14
                                                                      2
                                                rs540662886
                                                                            4979
15
                                                rs540662886
                                                                      1
                                                                             649
```

16				rs5406	662886		2 649
17	rs540662886			662886		1 3286	
18		rs540662886			2 3286		
	STRAND	FLAGS VA	RIANT_CLASS MINI	MISED SYMBO	OL_SOURCE HG	NC_I	D CANONICAL
1	1		SNV		HGNC	2402	23 YES
2	1		SNV		HGNC	2402	23 YES
3	-1		SNV		HGNC	2820)8
4	-1		SNV		HGNC	2820)8
5	1		SNV		HGNC	2528	34
6	1		SNV		HGNC	2528	34
7	1		SNV		HGNC	2528	34
8	1		SNV		HGNC	2528	34
9	1		SNV		HGNC	2528	34 YES
10	1		SNV		HGNC	2528	34 YES
11	-1		SNV		HGNC	2820	98 YES
12	-1		SNV			2820	
13	-1		SNV			2820	
14	-1		SNV		HGNC	2820)8
15	1		SNV			2528	
16	1		SNV			2528	
17		_start_NF	SNV			2528	
18		_start_NF	SNV			2528	
	TSL APPRIS			SWISSPROT	TREM		UNIPARC
1			ENSP00000343930	•	•		PI00001DFBF0
2		CCDS30550.1	ENSP00000343930	Q6TDP4	QOVGE6&B3KX		PI00001DFBF0
3			ENSP00000343864				JPI000022DAF4
4			ENSP00000343864				JPI000022DAF4
5			ENSP00000368717	Q494U1			PI00005764FF
6		CCDS53256.1	ENSP00000368717	Q494U1			PI00005764FF
7			ENSP00000368719	Q494U1			JPI0000D61E06
8			ENSP00000368719	Q494U1			JPI0000D61E06
9			ENSP00000368720	Q494U1			JPI00001416D8
10		CCDS4.1	ENSP00000368720	Q494U1	J3KS		PI00001416D8
11			ENSP00000414022	Q5SV97			JPI0000418FB0
12			ENSP00000414022	Q5SV97		Ţ	JPI0000418FB0
13							
14							
15							
16			THE PARAMETER AND THE				
17			ENSP00000462558				JPI000268AE1F
18			ENSP00000462558		J3KS	M5 U	JPI000268AE1F

```
GENE_PHENO SIFT PolyPhen DOMAINS HGVS_OFFSET
                                                          GMAF AFR_MAF AMR_MAF
1
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
2
3
4
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
5
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
6
                                                                             C:0
                                                     C:0.0008
                                                                    C:0
7
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
8
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
9
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
10
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
11
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
12
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
13
                                                                             C:0
                                                     C:0.0008
                                                                    C:0
14
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
15
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
16
                                                                             C:0
                                                     C:0.0008
                                                                    C:0
17
                                                                    C:0
                                                                             C:0
                                                     C:0.0008
18
                                                     C:0.0008
                                                                    C:0
                                                                             C:0
   EAS_MAF EUR_MAF SAS_MAF AA_MAF EA_MAF ExAC_MAF ExAC_Adj_MAF ExAC_AFR_MAF
1
       C:0 C:0.004
                         C:0
                                 C:0
                                                                       C:2.146e-04
                                                                  C:0
2
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
3
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
4
                                                                       C:2.146e-04
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
5
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
6
                                 C:0
                                                                       C:2.146e-04
       C:0 C:0.004
                         C:0
                                                                  C:0
       C:0 C:0.004
                                 C:0
                                                                       C:2.146e-04
                         C:0
                                                                  C:0
8
                                                                       C:2.146e-04
       C:0 C:0.004
                                 C:0
                                                                  C:0
                         C:0
9
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
10
       C:0 C:0.004
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
                         C:0
                                                                      C:2.146e-04
11
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
12
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
                                                                       C:2.146e-04
13
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
14
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
15
                                                                       C:2.146e-04
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
16
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
17
       C:0 C:0.004
                                                                       C:2.146e-04
                         C:0
                                 C:0
                                                                  C:0
18
       C:0 C:0.004
                         C:0
                                 C:0
                                                                  C:0
                                                                       C:2.146e-04
   ExaC_AMR_MAF ExaC_EAS_MAF ExaC_FIN_MAF ExaC_NFE_MAF ExaC_OTH_MAF
1
             C:0
                  C:0.0002281
                                  C:0.002986
                                                         C:0
                                                              C:1.606e-05
2
             C:0
                  C:0.0002281
                                  C:0.002986
                                                         C:0
                                                              C:1.606e-05
```

```
C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
4
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
5
                  C:0.0002281
                                 C:0.002986
                                                             C:1.606e-05
             C:0
                                                       C:0
6
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
7
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
8
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
9
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
10
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
             C:0
11
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
12
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
13
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
                                 C:0.002986
14
             C:0
                  C:0.0002281
                                                       C:0
                                                             C:1.606e-05
15
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
16
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
17
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
18
             C:0
                  C:0.0002281
                                 C:0.002986
                                                       C:0
                                                             C:1.606e-05
   ExAC_SAS_MAF CLIN_SIG SOMATIC PHENO PUBMED MOTIF_NAME MOTIF_POS HIGH_INF_POS
             C:0
1
2
             C:0
3
             C:0
4
             C:0
5
             C:0
6
             C:0
7
             C:0
8
             C:0
9
             C:0
10
             C:0
11
             C:0
12
             C:0
13
             C:0
14
             C:0
15
             C:0
16
             C:0
17
             C:0
18
             C:0
   MOTIF_SCORE_CHANGE LoF LoF_filter LoF_flags LoF_info"
1
2
3
4
```

```
6
7
8
9
10
11
12
13
14
15
16
17
18
>
```

2.12 Print a SNPEFF table for a Variant

Code:

```
# load rbcf
library(rbcf)
# A vcf
filename <- "../tests/data/rotavirus_rf.ann.vcf.gz"
# we don't need the index for this file
fp <- bcf.open(filename,FALSE)</pre>
# current variant
vc <- NULL
while(!is.null(vc<-bcf.next(fp))) {</pre>
        #find the first variant having an INFO/ANN attribute
        if(variant.has.attribute(vc,"ANN")) break;
if(!is.null(vc)) {
        # get SNPEFF table
        predictions<-variant.snpeff(vc)</pre>
# dispose the vcf reader
bcf.close(fp)
# show
predictions
```

```
> # load rbcf
> library(rbcf)
> # A vcf
> filename <- "../tests/data/rotavirus_rf.ann.vcf.gz"</pre>
> # we don't need the index for this file
> fp <- bcf.open(filename,FALSE)</pre>
> # current variant
> vc <- NULL
> while(!is.null(vc<-bcf.next(fp))) {</pre>
        #find the first variant having an INFO/ANN attribute
        if(variant.has.attribute(vc,"ANN")) break;
        }
> if(!is.null(vc)) {
        # get SNPEFF table
       predictions<-variant.snpeff(vc)</pre>
> # dispose the vcf reader
> bcf.close(fp)
[1] TRUE
> # show
> predictions
               Annotation Annotation_Impact
  Allele
                                                Gene_Name
       C missense_variant
                                    MODERATE Gene 18 3284 Gene 18 3284
 Feature_Type Feature_ID Transcript_BioType Rank
                                                     HGVS.c
1 transcript AAA47319.1
                            protein_coding 1/1 c.952A>C p.Lys318Gln
  cDNA.pos / cDNA.length CDS.pos / CDS.length AA.pos / AA.length Distance
                952/3267
                                                         318/1088
                                      952/3267
 ERRORS / WARNINGS / INFO'"
```